

Catheter Spigot - AECS9004

Product Data Sheet

Description	Sterile, single use, closures for catheters or tubes.
Intended Use	A sterile plug-like device intended to be inserted into a catheter to prevent flow through. It may also be referred to as a catheter spigot or catheter tip occluder. (GMDN 31667)
Duration of Use	"Short Term" : Normally intended for continuous use of more than 60 minutes but less than 30 days.
Materials	PP - Polypropylene
Product Information	Tip OD: 5.6 mm Length: 54.5 mm Colour: Green Universal size is moulded to fit all sizes of catheters or tubes Graduated design easily connects tubing of various internal diameters Sterile, individual peel pouches Latex Free
Packaging	500 units in the carton
Quality System	ISO 13485
EU Classification	Class Is
Sterilised by	Ethylene Oxide
Shelf Life	5 years
Storage & Safety	Toxicity: Manufactured from biocompatible materials. Storage: Store in a clean, cool, dry place away from direct sunlight.



STERILE EO



GCMEDICA ENTERPRISE LTD., (WUXI)
Loujin Industrial Park, Shuofang, Wuxi 214143, Jiangsu, China



Kingsmead Service B.V.
Zonnehof 36, 2632 BE, Nootdorp, Netherland

© 2025 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.