



Product Information



Sterile hypodermic needle for single use

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect¹.

Description	Raw materials
1 Hub	Polypropylene + master batch
2 Cannula	Stainless steel
Adhesive	Epoxy glue
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.



Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶.



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the packaging
- Do not store at extreme temperature and humidity



Precautions

- Check integrity of the package before use
- Do not reuse
- Do not use if unit package is damaged
- Dispose of safely after single use to avoid risk of infection

¹ Data on file

Your needles of choice for a wide range of applications

Product range

Product Code	Colour	Cannula ø	Cannula Length	Bevel
AN*1838R1	Pink	18G 1.2 mm	1 ½ inch 38 mm	Regular 11°
AN*1838S1	Pink	18G 1.2 mm	1 ½ inch 38 mm	Short 17°
AN*1925R1	Cream	19G 1.1 mm	1 inch 25 mm	Regular 11°
AN*1938R1	Cream	19G 1.1 mm	1 ½ inch 38 mm	Regular 11°
AN*1938S1	Cream	19G 1.1 mm	1 ½ inch 38 mm	Short 17°
AN*1950R1	Cream	19G 1.1 mm	2 inch 50 mm	Regular 11°
AN*2025R1	Yellow	20G 0.9 mm	1 inch 25 mm	Regular 11°
AN*2038R1	Yellow	20G 0.9 mm	1 ½ inch 38 mm	Regular 11°
AN*2116R1	Green	21G 0.8 mm	⅝ inch 16 mm	Regular 11°
AN*2125R1	Green	21G 0.8 mm	1 inch 25 mm	Regular 11°
AN*2138R1	Green	21G 0.8 mm	1 ½ inch 38 mm	Regular 11°
AN*2150R1	Green	21G 0.8 mm	2 inch 50 mm	Regular 11°
AN*2232R1	Black	22G 0.7 mm	1 ¼ inch 32 mm	Regular 11°
AN*2238R1	Black	22G 0.7 mm	1 ½ inch 38 mm	Regular 11°
AN*2316R1	Blue	23G 0.6 mm	⅝ inch 16 mm	Regular 11°
AN*2325R1	Blue	23G 0.6 mm	1 inch 25 mm	Regular 11°
AN*2332R1	Blue	23G 0.6 mm	1 ¼ inch 32 mm	Regular 11°
AN*2425R1	Purple	24G 0.55 mm	1 inch 25 mm	Regular 11°
AN*2516R1	Orange	25G 0.5 mm	⅝ inch 16 mm	Regular 11°
AN*2525R1	Orange	25G 0.5 mm	1 inch 25 mm	Regular 11°
AN*2613R1	Brown	26G 0.45 mm	½ inch 13 mm	Regular 11°
AN*2623R1	Brown	26G 0.45 mm	⅞ inch 23 mm	Regular 11°
AN*2716R1	Grey	27G 0.4 mm	⅝ inch 16 mm	Regular 11°
AN*2719R1	Grey	27G 0.4 mm	¾ inch 19 mm	Regular 11°
AN*3013R1	Yellow	30G 0.3 mm	½ inch 13 mm	Regular 11°

Packaging

Single unit Blister pack (film/paper)

Unit box 100 pieces

Dimensions 105 x 84 x 86 mm

For item numbers AN*1950R1 and

AN*2150R1 dimensions 105 x 84 x 110 mm

Shipping carton 10.000 pieces

Dimensions 450 x 440 x 450 mm

For item numbers AN*1950R1 and

AN*2150R1 5.000 pieces,

dimensions 545 x 240 x 450 mm

Regulatory Status

Class IIa, sterile product following

MDD 93/42/EEC

Notified body TÜV Product Service

GmbH, TÜV SUD Group 0123

Manufacturer Zhejiang Kindly Medical

Devices Co. Ltd. (China)

EC representative Shanghai International

Corp. GmbH (Europe) - Hamburg, Germany

Distributor TERUMO EUROPE N.V.

Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a wide range from 18G up to 30G

The outside diameter of the needle can be recognized by the ISO colour-coded hub.



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Your Terumo partner